

Nur für Forschungszwecke

AgBR1 Monoklonaler Antikörper

Katalog-Nr.: 67743-1-Ig



Allgemeine Informationen

Katalog-Nr.: 67743-1-Ig	GenBank-Zugangsnummer: XM_001844620	Reinigungsmethode: Protein-G-Reinigung
Größe: 150ul, Konzentration: 1000 µg/ml von 6034353	GeneID (NCBI): 6034353	CloneNo.: 2D1F10
Wirt: Maus	Vollständiger Name: AgBR1	Empfohlene Verdünnungen: WB 1:20000-1:100000
Isotyp: IgG1	Berechnete Masse: 49 kDa	
Immunogen Katalognummer: AG30168	Beobachtete Masse: 40 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : Ag30168,
Getestete Reaktivität: Stechmücke	

Hintergrundinformationen

AgBR1, a mosquito salivary protein also known as A. aegypti bacteria-responsive protein 1, is found to induce inflammatory responses at the bite site (PMID: 32218189). Various studies suggest that AgBR1 plays a role in transmitting or facilitating infections of several viruses via Aedes aegypti mosquitoes (PMID: 30858571, 31312526).

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 67743-1-Ig at dilution of 1:160000 incubated at room temperature for 1.5 hours.