

Nur für Forschungszwecke

INADL Monoklonaler Antikörper

Katalog-Nr.: 67833-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67833-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml von 10207

Nanodrop;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG15953

GenBank-Zugangsnummer:

BC142661

GeneID (NCBI):

von 10207

Vollständiger Name:

InaD-like (Drosophila)

Berechnete Masse:

1801 aa, 196 kDa

Beobachtete Masse:

196 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

3C6A4

Empfohlene Verdünnungen:

WB 1:2000-1:10000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human, Ratte

Positivkontrollen:

WB: LNCaP-Zellen, HEK-293-Zellen, HeLa-Zellen,
Jurkat-Zellen, PC-12-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 67833-1-Ig (INADL antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.