

Nur für Forschungszwecke

GST Tag Rekombinanter Antikörper

Katalog-Nr.:80006-1-RR



Allgemeine Informationen

Katalog-Nr.: 80006-1-RR	GenBank-Zugangsnummer: FN315687	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von Nanodrop;	GeneID (NCBI): Vollständiger Name:	CloneNo.: 9B24
Wirt: Kaninchen	Berechnete Masse: 26 kDa	Empfohlene Verdünnungen: WB 1:5000-1:50000
Isotyp: IgG		
Immunogen Katalognummer: AG0040		

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen: WB : Recombinant protein,
Getestete Reaktivität: Rekombinantes Protein	

Hintergrundinformationen

GST (Glutathione S-Transferase) is a widely used protein tag encoded by the Schistosoma japonicum. GST provides both an easily detectable tag and a simple purification process with little effect on the biological function of the protein of interest. Antibodies to GST are useful for checking protein expression both in plaques and on western blots as well as for immunoaffinity purification of proteins expressed in cells. This antibody is a Rabbit recombinant antibody to GST tag.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



GST tagged protein and other negative samples were subjected to SDS PAGE followed by western blot with 80006-1-RR (GST Tag antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.