

Nur für Forschungszwecke

TTYH1 Rekombinanter Antikörper

Katalog-Nr.:80263-1-RR



Allgemeine Informationen

Katalog-Nr.: 80263-1-RR	GenBank-Zugangsnummer: BC011347	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von57348	GeneID (NCBI): 57348	CloneNo.: 8G21
Nanodrop;	Vollständiger Name: tweety homolog 1 (Drosophila)	Empfohlene Verdünnungen: WB 1:5000-1:50000
Wirt: Kaninchen	Berechnete Masse: 49-50 kDa	
Isotyp: IgG		
Immunogen Katalognummer: AG25290		

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen:
Getestete Reaktivität: Hausschwein, Human, Kaninchen	WB : Hausschwein-Cerebellum-Gewebe, humanes Hodengewebe, Kaninchen-Cerebellum-Gewebe, Kaninchenhirngewebe

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 80263-1-RR (TTYH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.