

Nur für Forschungszwecke

FGFR4 Rekombinanter Antikörper

Katalog-Nr.: 81070-1-RR



Allgemeine Informationen

Katalog-Nr.: 81070-1-RR	GenBank-Zugangsnummer: BC011847	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul, Konzentration: 2292 µg/ml von 2264 Nanodrop;	GeneID (NCBI): 2264	CloneNo.: 2M22
Wirt: Kaninchen	Vollständiger Name: fibroblast growth factor receptor 4	Empfohlene Verdünnungen: WB 1:4000-1:40000
Isotyp: IgG	Berechnete Masse: 88 kDa	
Immunogen Katalognummer: AG31013	Beobachtete Masse: 100-110 kDa	

Anwendungen

Geprüfte Anwendungen: WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human	WB : MOLT-4-Zellen, A549-Zellen, HeLa-Zellen, HepG2-Zellen, HT-29-Zellen, Jurkat-Zellen, Raji-Zellen

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 81070-1-RR (FGFR4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.