

Allgemeine Informationen

Katalog-Nr.: 81471-1-RR	GenBank-Zugangsnummer: NM_001080432	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul, Konzentration: 1000 µg/ml von 79068	GeneID (NCBI): von 79068	CloneNo.: 4L4
Nanodrop;	Vollständiger Name: fat mass and obesity associated	Empfohlene Verdünnungen: WB 1:2000-1:10000 IHC 1:500-1:2000
Wirt: Kaninchen	Berechnete Masse: 58 kDa	
Isotyp: IgG	Beobachtete Masse: 58 kDa	
Immunogen Katalognummer: AG26095		

Anwendungen

Geprüfte Anwendungen: IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human	WB : HeLa-Zellen, HEK-293T-Zellen, Jurkat-Zellen, LNCaP-Zellen
Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.	IHC : humanes Leberkarzinomgewebe,


Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 81471-1-RR (FTO antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using 81471-1-RR (FTO antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).