

Nur für Forschungszwecke

# SARS-CoV-2 Spike Rekombinanter Antikörper

Katalog-Nr.: 91351-PTG



## Allgemeine Informationen

Katalog-Nr.:  
91351-PTG

Größe:  
100ug

Wirt:  
Mensch

Isotyp:  
IgG1

GenBank-Zugangsnummer:  
NC\_045512

GeneID (NCBI):  
43740568

Vollständiger Name:  
SARS-CoV-2 Spike Protein

Berechnete Masse:  
141 kDa

Reinigungsmethode:  
Protein-A-Chromatographie

CloneNo.:  
AM006415

Empfohlene Verdünnungen:  
Sample dependent. To be determined by the end user.

## Anwendungen

Geprüfte Anwendungen:  
WB, ELISA

Getestete Reaktivität:  
Virus

## Hintergrundinformationen

COVID-19, which is short for coronavirus disease 2019, is the official name of the respiratory disease caused by infection with the novel coronavirus SARS-CoV-2. The virus that causes COVID-19 was named SARS-CoV-2 because it is a coronavirus genetically similar to, yet distinct from, the virus that caused the severe acute respiratory syndrome (SARS) outbreak in 2003. Studying the details of how this virus replicates and causes the disease will allow scientists and physicians to more rapidly develop fast and accurate methods of detection as well as to deploy therapeutic and vaccine strategies. This antibody was derived from COVID-19 patients who have cleared the virus. Patient serum IgG was sequenced and expressed as full-length IgG1 with human immunoglobulin heavy and light chains in mammalian 293 cells.

## Lagerung

Lagerungsbedingungen:

Avoid repeated freeze/thaw cycles by aliquoting items into single-use fractions for storage at -20°C for up to 2 years. Keep all reagents on ice when not in storage. This product is guaranteed for 12 months from date of receipt.

Lagerungspuffer:

140 mM HEPES, pH 7,5, 70 mM NaCl, 32 mM NaOAc, 0,035% Natriumazid und 30% Glycerin.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

\*\*\* 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is made by Active Motif® and available through Proteintech Group, Inc.

## Ausgewählte Validierungsdaten



SARS-CoV-2 Spike Antibody (AM006415) tested by Western blot. Lysates from HEK293 cells transfected (Lane 1 and 3) or untransfected (Lane 2) were probed with 0.5 µg/mL of rAb (Lane 1 & 2) or an anti-flag Ab (Lane 3). Lane 1: transfected flag-tagged SARS-CoV-2 S1 subunit. Lane 2: untransfected. Lane 3: transfected flag-tagged SARS-CoV-2 S1 subunit. Data provided by Active Motif®.



SARS-CoV-2 Spike Antibody (clone AM006415) tested by ELISA. SARS-CoV-2 Spike RBD protein was coated onto microtiter plates at 0.5 µg/mL and then incubated with a dilution series of SARS-CoV-2 Spike Antibody (clone AM009105). Bound antibodies were detected with anti-human IgG conjugated to horseradish peroxidase (HRP) followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm. Data provided by Active Motif®.