

Nur für Forschungszwecke

# CD4 Monoklonaler Antikörper

Katalog-Nr.: AB-65134



## Allgemeine Informationen

<b>Katalog-Nr.:</b> AB-65134	<b>GenBank-Zugangsnummer:</b> BC025782	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests, 5 µl/test	<b>GeneID (NCBI):</b> 920	<b>CloneNo.:</b> OKT4
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> CD4 molecule	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 404 nm / 458 nm
<b>Isotyp:</b> IgG2b, kappa	<b>Berechnete Masse:</b> 55 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

## Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Lagerungspuffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



$1 \times 10^6$  human PBMCs were surface co-stained with APC Anti-Human CD3 and 5  $\mu$ l Atlantic Blue™ Anti-Human CD4 (AB-65134, Clone: OKT4) or Mouse IgG2b Isotype Control. Cells were not fixed. Lymphocytes were gated.



$1 \times 10^6$  human PBMCs were surface stained with 5  $\mu$ l Atlantic Blue™ Anti-Human CD4 (AB-65134, Clone: OKT4) or Mouse IgG2b Isotype Control. Cells were not fixed. Lymphocytes were gated.