

Nur für Forschungszwecke

# CD29 Monoklonaler Antikörper

Katalog-Nr.: AB-65191



## Allgemeine Informationen

<b>Katalog-Nr.:</b> AB-65191	<b>GenBank-Zugangsnummer:</b> BC020057	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100tests, 5 µl/test	<b>GeneID (NCBI):</b> 3688	<b>CloneNo.:</b> TS2/16
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> integrin, beta 1 (fibronectin receptor, beta polypeptide, antigen CD29 includes MDF2, MSK12)	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 404 nm / 458 nm
<b>Isotyp:</b> IgG1, kappa	<b>Berechnete Masse:</b> 88 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

Integrin beta-1 (ITGB1), also named as CD29, is a 130 kDa single chain type I glycoprotein that is expressed in a heterodimeric complex with one of six distinct  $\alpha$  subunits, comprising the very late activation antigen (VLA) subfamily of adhesion receptors. It is one of the essential surface molecules expressed on human MSC from bone marrow and other sources. The  $\beta$ 1 subunit is also broadly expressed on lymphocytes and monocytes, weakly expressed on granulocytes, and not expressed on erythrocytes. These receptors are involved in a variety of cell-cell and cell-matrix interactions.

## Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Lagerungspuffer:**  
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD29 (AB-65191, Clone:TS2/16) or unstained. Cells were not fixed. Monocytes were gated.



1X10<sup>6</sup> human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD29 (AB-65191, Clone:TS2/16) or unstained. Cells were not fixed. Monocytes were gated.