

Nur für Forschungszwecke

CD71 Monoklonaler Antikörper

Katalog-Nr.: AB-65235



Allgemeine Informationen

Katalog-Nr.:

AB-65235

Größe:

100tests, 5 µl/test

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC001188

GeneID (NCBI):

7037

Vollständiger Name:

transferrin receptor (p90, CD71)

Berechnete Masse:

85 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

OKT9

Anregungs-/Emissionsmaxima-
Wellenlängen:

404 nm / 458 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ PHA-activated human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD71 (AB-65235, Clone:OKT9) (red) or Mouse IgG1 Isotype Control (blue). Unactivated human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD71 (black). Cells were not fixed. Lymphocytes were gated.