

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: AB-65244

Allgemeine Informationen

Katalog-Nr.:

AB-65244

Größe:

100µg, 200 µg/ml

Wirt:

Maus

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

MOPC-173

Anregungs-/Emissionsmaxima-Wellenlängen:

404 nm / 458 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 0.2 ug Atlantic Blue™ Anti-Mouse IgG2a Isotype Control (AB-65244, Clone:MOPC-173) (red), Atlantic Blue™ Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.