

Nur für Forschungszwecke

CD35 Monoklonaler Antikörper

Katalog-Nr.: AB-65254



Allgemeine Informationen

Katalog-Nr.:
AB-65254

Größe:
100tests, 5 µl/test

Wirt:
Maus

Isotyp:
IgG1, kappa

GenBank-Zugangsnummer:
NM_000573

GeneID (NCBI):
1378

Vollständiger Name:
complement component (3b/4b)
receptor 1 (Knops blood group)

Berechnete Masse:
224 kDa

Reinigungsmethode:
Affinitätsreinigung

CloneNo.:
E11

Anregungs-/Emissionsmaxima-
Wellenlängen:
404 nm / 458 nm

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD35 (AB-65254, Clone:E11) or 5 ul Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 5 ul Atlantic Blue™ Anti-Human CD35 (AB-65254, Clone:E11) or 5 ul Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.