

Nur für Forschungszwecke

CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.: APC-65095



Allgemeine Informationen

Katalog-Nr.: APC-65095	GenBank-Zugangsnummer: BC036343	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug, 0.2 mg/ml	GeneID (NCBI): 12978	CloneNo.: AFS98
Wirt: Ratte	Vollständiger Name: colony stimulating factor 1 receptor	Anregungs-/Emissionsmaxima-Wellenlängen: 650 nm / 660 nm
Isotyp: IgG2a, kappa		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 0.13 ug APC Anti-Mouse CD115 (c-fms) (APC-65095, Clone: AFS98) (blue) or 0.13 ug isotype control (black). Cells were not fixed.