

Nur für Forschungszwecke

CD24 Monoklonaler Antikörper

Katalog-Nr.: APC-65121



Allgemeine Informationen

Katalog-Nr.: APC-65121	GenBank-Zugangsnummer: BC075622	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug, 0.2 mg/ml	GeneID (NCBI): 12484	CloneNo.: M1/69
Wirt: Ratte	Vollständiger Name: CD24a antigen	Anregungs-/Emissionsmaxima- Wellenlängen: 650 nm / 660 nm
Isotyp: IgG2b		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.
Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with FITC Anti-Mouse CD3 (FITC-65077, Clone: 17A2) and 0.06 ug APC-Rat IgG2b isotype control (left) or 0.06 ug APC Anti-Mouse CD24 (APC-65121, Clone: M1/69) (right). Cells were not fixed.