

Nur für Forschungszwecke

CD11c Monoklonaler Antikörper

Katalog-Nr.: APC-65130

3 Publikationen



Allgemeine Informationen

Katalog-Nr.:

APC-65130

Größe:

100ug, 0.2 mg/ml

Wirt:

Armenian Hamster

Isotyp:

IgG

GenBank-Zugangsnummer:

BC167225

GeneID (NCBI):

16411

Vollständiger Name:

integrin alpha X

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

N418

Anregungs-/Emissionsmaxima-
Wellenlängen:

650 nm / 660 nm

Anwendungen

Geprüfte Anwendungen:

FC

In Publikationen genannte Anwendungen:

FC

Getestete Reaktivität:

Maus

Zitierte Arten:

Maus

Hintergrundinformationen

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits (PMID: 9779984). CD11c, also known as integrin αX , is a type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, a subset of B cells, NK cells and dendritic cells (PMID: 2897326; 1680915; 1694698; 17389580). As a result of its high level of expression on most dendritic cells, CD11c is typically considered to be a marker of conventional dendritic cells (PMID: 27119555). CD11c forms an α/β heterodimer with CD18 (integrin $\beta 2$). CD11c/CD18 acts a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis (PMID: 1671533).

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Xiaojin Chen	36297436	Pharmaceutics	FC
Dingyi Wang	35236203	Drug Deliv	FC
Xudong Fan	36382303	Asian J Pharm Sci	FC

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.1% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1x10⁶ mouse splenocytes were surface co-stained with FITC Anti-Mouse MHC Class II (I-A/I-E) and 0.2 ug APC Anti-Mouse CD11c (APC-65130, Clone: N418) or 0.2 ug Isotype Control. Cells were not fixed.