

Nur für Forschungszwecke

CD1a Monoklonaler Antikörper

Katalog-Nr.: APC-65172



Allgemeine Informationen

Katalog-Nr.:

APC-65172

Größe:

100tests, 10 µl/test

Wirt:

Maus

Isotyp:

Mouse IgG1

GenBank-Zugangsnummer:

BC031645

GeneID (NCBI):

909

Vollständiger Name:

CD1a molecule

Berechnete Masse:

327 aa, 37 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit Allophycocyanin (APC) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:

HI149

Anregungs-/Emissionsmaxima-Wellenlängen:

650 nm / 660 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ MOLT-4 cells were surface stained with 10.00 ul/test APC-Anti-Human CD1a antibody (APC-65172, Clone: HI 149) (red) or APC-mouse IgG1 isotype control (blue). Cells were not fixed.