

Nur für Forschungszwecke

CD68 Monoklonaler Antikörper

Katalog-Nr.: APC-65187



Allgemeine Informationen

Katalog-Nr.:

APC-65187

Größe:

100tests, 10 µl/test

Wirt:

Maus

Isotyp:

Mouse IgG2a

GenBank-Zugangsnummer:

BC015557

GeneID (NCBI):

968

Vollständiger Name:

CD68 molecule

Berechnete Masse:

37 kDa

Reinigungsmethode:

Der gereinigte Antikörper ist mit Allophycocyanin (APC) unter optimalen Bedingungen konjugiert. Das Konjugat wird mit Größenausschluss-Chromatographie gereinigt.

CloneNo.:

Y1/82A

Anregungs-/Emissionsmaxima-Wellenlängen:

650 nm / 660 nm

Anwendungen

Geprüfte Anwendungen:

FC (Intra)

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood monocytes were intracellularly stained with 10.00 uL/test APC Anti-Human CD68 (APC-65187, Clone:Y1/82A) (red) or isotype control antibody (blue). Cells were fixed with 4% PFA and permeabilized with 0.1% Triton X-100.