

Nur für Forschungszwecke

# IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: APC-65209

## Allgemeine Informationen

Katalog-Nr.:

APC-65209

Größe:

100ug, 0.2 mg/mL

Wirt:

Ratte

Isotyp:

IgG2a, κ

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

2A3

Anregungs-/Emissionsmaxima-Wellenlängen:

650 nm / 660 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

## Hintergrundinformationen

The 2A3 immunoglobulin is useful as an isotype-matched control. The 2A3 immunoglobulin binds to TNP (trinitrophenol), a hapten not expressed on mouse or human cells. It is used as an isotype control for rat IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

## Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

Phosphatbasierter Puffer mit 0,09% Natriumazid und 0,1% Gelatine, pH 7,2.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)

W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> BALB/c mouse splenocytes were surface stained with 0.5 ug APC Rat IgG2a Isotype Control (APC-65209, Clone: 2A3) (blue) or 0.5 ug APC Anti-Mouse CD4 (APC-65141, Clone: RM4-5) (red). Cells were not fixed.