

RSPO2 Polyklonaler Antikörper

Katalog-Nr.: Biotin-17781

Allgemeine Informationen

Katalog-Nr.:

Biotin-17781

Größe:

100ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Kaninchen

Isotyp:

IgG

Immunogen Katalognummer:

AG11216

GenBank-Zugangsnummer:

BC027938

GeneID (NCBI):

340419

Vollständiger Name:

R-spondin 2 homolog (Xenopus
laevis)

Berechnete Masse:

243 aa, 28 kDa

Beobachtete Masse:

33 kDa

Reinigungsmethode:

Antigen-Affinitätsreinigung

Anwendungen

Geprüfte Anwendungen:

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

RSPO2 is one of the R-spondins (RSPOs) which comprise a family of four secreted proteins that have important roles in cell proliferation, cell fate determination and organogenesis. RSPOs typically exert their effects by potentiating the Wnt/beta-catenin signaling pathway. Overexpression of Rspo2 is sufficient to promote bone morphogenetic protein-induced mineralization Rspo2 was identified as a novel common integration site (CIS) for the mouse mammary tumor virus (MMTV) in viral induced mouse mammary tumors.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

