

Nur für Forschungszwecke

CD45 Monoklonaler Antikörper

Katalog-Nr.: CL405-65228



Allgemeine Informationen

Katalog-Nr.:

CL405-65228

Größe:

100ug, 0.5 mg/ml

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

XM_006249910

GeneID (NCBI):

24699

Vollständiger Name:

protein tyrosine phosphatase,
receptor type, C

Berechnete Masse:

143 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

OX1

Anregungs-/Emissionsmaxima-
Wellenlängen:

399 nm / 422 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Ratte

Hintergrundinformationen

CD45, also known as protein tyrosine phosphatase, receptor type C, is a type I transmembrane protein expressed on the surface of all haematopoietic cells with the exception of erythrocytes and platelets (PMID: 3489673; 28615666). CD45 is a pan-haematopoietic cell marker and has been shown to be essential for T- and B-cell activation and signalling (PMID: 9429890; 16378097).

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Rat CD45 (CL405-65228, Clone: OX1) (red), or 0.5 ug Control Antibody. Cells were not fixed.