

Nur für Forschungszwecke

CD71 Monoklonaler Antikörper

Katalog-Nr.:CL405-65235



Allgemeine Informationen

Katalog-Nr.: CL405-65235	GenBank-Zugangsnummer: BC001188	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 7037	CloneNo.: OKT9
Wirt: Maus	Vollständiger Name: transferrin receptor (p90, CD71)	Anregungs-/Emissionsmaxima-Wellenlängen: 399 nm / 422 nm
Isotyp: IgG1, kappa	Berechnete Masse: 85 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ PHA-treated (day 3) human PBMCs were surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD71 (CL405-65235, Clone:OKT9) (red) or unstained (blue). The black dashed line indicates unactivated PBMCs surface stained with 5 ul CoraLite® Plus 405 Anti-Human CD71. Cells were not fixed. Lymphocytes were gated.