

Nur für Forschungszwecke

CD43 Monoklonaler Antikörper

Katalog-Nr.:CL405-65237



Allgemeine Informationen

Katalog-Nr.:
CL405-65237

Größe:
100ug, 500 µg/ml

Wirt:
Maus

Isotyp:
IgG1, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):
24796

Vollständiger Name:
sialophorin

Berechnete Masse:
38 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

W3/13

Anregungs-/Emissionsmaxima-
Wellenlängen:

399 nm / 422 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Wistar rat splenocytes were surface stained with 0.5 ug CoraLite® Plus 405 Anti-Rat CD43 (CL405-65237, Clone: W3/13) or 0.5 ug Control Antibody. Cells were not fixed.