

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL405-65244

Allgemeine Informationen

Katalog-Nr.:
CL405-65244
Größe:
100ug , 200 µg/ml
Wirt:
Maus
Isotyp:
IgG2a, kappa

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
Affinitätsreinigung
CloneNo.:
MOPC-173
Anregungs-/Emissionsmaxima-Wellenlängen:
399 nm / 422 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 human PBMCs were surface stained with 0.2 μ g Coralite® Plus 405 Anti-Mouse IgG2a Isotype Control (CL405-65244, Clone: MOPC-173) (red), CL405 Anti-Human CD5 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.