

Nur für Forschungszwecke

Mitofilin Polyklonaler Antikörper

Katalog-Nr.:CL488-10179

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.: CL488-10179	GenBank-Zugangsnummer: BC002412	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul , Konzentration: 1000 µg/ml von10989 Nanodrop;	GeneID (NCBI): Vollständiger Name: inner membrane protein, mitochondrial (mitofilin)	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Wirt: Kaninchen	Berechnete Masse: 90 kDa	
Isotyp: IgG	Beobachtete Masse: 80-90 kDa	
Immunogen Katalognummer: AG0102		

Anwendungen

Geprüfte Anwendungen:
FC (Intra)

Getestete Reaktivität:
Human, Maus, Ratte

Hintergrundinformationen

IMMT (also known as mitofilin) is an inner mitochondrial membrane protein that is preferentially expressed in heart tissue and is commonly used as the marker for mitochondria. Mitofilin is known to be a critical organizer of mitochondrial cristae morphology and is indispensable for normal mitochondrial function. Three isoforms of mitofilin exist due to the alternative splicing. All three proteins of 87kD, 89kD and 80kD can be detected in immunoblot analysis with this antibody.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ HEK-293T cells were intracellularly stained with 0.4 ug CoraLite® Plus 488 Anti-Human Mitofilin (CL488-10179) (red), or 0.4 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).