

Nur für Forschungszwecke

C1QBP Polyklonaler Antikörper

Katalog-Nr.:CL488-24474

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.: CL488-24474	GenBank-Zugangsnummer: BC013731	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul , Konzentration: 1000 µg/ml von708 Nanodrop;	GeneID (NCBI): Vollständiger Name: complement component 1, q subcomponent binding protein	Empfohlene Verdünnungen: IF 1:50-1:500
Wirt: Kaninchen	Berechnete Masse: 282 aa, 31 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Isotyp: IgG	Beobachtete Masse: 32 kDa	
Immunogen Katalognummer: AG19773		

Anwendungen

Geprüfte Anwendungen: FC (Intra), IF	Positivkontrollen: IF : NIH/3T3-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

C1QBP, also named as gC1q receptor (gC1qR), p32, p33, and hyaluronan-binding protein 1 (HABP1), is a protein initially copurified with splicing factor SF2 (PMID: 1830244). The protein is synthesized as a pro-protein of 282 amino acids (aa) that is post-translationally processed by removal of the initial 73 aa to a mature protein of 209 aa (PMID: 8262387). C1QBP is an evolutionary conserved and ubiquitously expressed multifunctional protein and has been reported to be a predominantly mitochondrial matrix protein involved in inflammation and infection processes, mitochondrial ribosome biogenesis, regulation of apoptosis and nuclear transcription, and pre-mRNA splicing (PMID: 28942965).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed NIH/3T3 cells using CoraLite® Plus 488 C1QBP antibody (CL488-24474) at dilution of 1:100.



1×10^6 HeLa cells were intracellularly stained with 0.4 μg CoraLite® Plus 488 Anti-Human C1QBP (CL488-24474) (red), or 0.4 μg Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).