

Nur für Forschungszwecke

CD11c/Integrin Alpha X Monoklonaler Antikörper



Katalog-Nr.:CL488-60258

1 Publikationen

Allgemeine Informationen

Katalog-Nr.: CL488-60258	GenBank-Zugangsnummer: BC038237	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von3687 Nanodrop;	GeneID (NCBI): von3687	CloneNo.: 2F1C10
Wirt: Maus	Vollständiger Name: integrin, alpha X (complement component 3 receptor 4 subunit)	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG2a	Berechnete Masse: 1169 aa, 129 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Immunogen Katalognummer: AG11350		

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : humanes Tonsillitisgewebe,
In Publikationen genannte Anwendungen: IF	
Getestete Reaktivität: Human	
Zitierte Arten: Maus	

Hintergrundinformationen

Integrins are cell adhesion receptors that are heterodimers composed of non-covalently associated α and β subunits. CD11c/Integrin alpha X is a 145-150 kDa type I transmembrane glycoprotein present on a variety of cells, including monocytes/macrophages, granulocytes, NK cells and dendritic cells. Integrin alpha X/beta 2 acts a receptor for fibrinogen and is important in monocyte adhesion and chemotaxis. And this antibody is CL488(Ex/Em 488 nm/515 nm) conjugated.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Kirsten M Reeves	34615724	Clin Cancer Res	IF

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 CD11c/Integrin Alpha X antibody (CL488-60258, Clone: 2F1C10) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 CD11c/Integrin Alpha X antibody (CL488-60258, Clone: 2F1C10) at dilution of 1:200.