

Nur für Forschungszwecke

CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.:CL488-65095



Allgemeine Informationen

Katalog-Nr.:

CL488-65095

Größe:

100tests , 5 µl/test

Wirt:

Ratte

Isotyp:

Rat IgG2a, kappa

GenBank-Zugangsnummer:

BC036343

GeneID (NCBI):

12978

Vollständiger Name:

colony stimulating factor 1 receptor

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

AFS98

Anregungs-/Emissionsmaxima-

Wellenlängen:

493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern. Vor Licht schützen.

Lagerungspuffer:

PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ Balb/c mouse peritoneal macrophages were surface stained with 5 ul CoraLite® Plus 488-conjugated Anti-Mouse CD115 (c-fms) (CL488-65095, Clone: AF598) (green) or CoraLite® Plus 488-conjugated Rat IgG2a isotype control (black). Cells were not fixed.