

Nur für Forschungszwecke

CD19 Monoklonaler Antikörper

Katalog-Nr.:CL488-65110



Allgemeine Informationen

| | | |
|---------------------------------------|---|--|
| Katalog-Nr.: CL488-65110 | GenBank-Zugangsnummer: BC006338 | Reinigungsmethode: Affinitätsreinigung |
| Größe: 100tests , 5 µl/test | GeneID (NCBI): 930 | CloneNo.: H1B19 |
| Wirt: Maus | Vollständiger Name: CD19 molecule | Anregungs-/Emissionsmaxima-Wellenlängen: 493 nm / 522 nm |
| Isotyp: IgG1, kappa | Berechnete Masse: 556 aa, 61 kDa | |

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

CD19 is a 95 kDa type I transmembrane glycoprotein belonging to the immunoglobulin superfamily (PMID: 2472450). It is expressed by B cells and follicular dendritic cells. CD19 is up-regulated at the step of B-lineage commitment during the differentiation of the hematopoietic stem cell, it remains on during subsequent stages of differentiation until finally down-regulated during terminal differentiation into plasma cells (PMID: 8528044). CD19 is involved in B cell development, activation and differentiation. It is the dominant component for the signaling complex on B cells that includes CD21 (CR2), CD81 (TAPA-1) and CD225 and acts as a critical co-receptor for BCR signal transduction (PMID: 23210908).

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood lymphocytes were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Human CD3 (CL647-65151, Clone: UCHT1) and 5 ul CoraLite® Plus 488-conjugated Anti-Human CD19 (CL488-65110, Clone: HIB19). Cells were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with 5 ul CoraLite® Plus 647-conjugated Anti-Human CD3 (CL647-65151, Clone: UCHT1) and CoraLite® Plus 488-conjugated Mouse IgG1 isotype control. Cells were not fixed.