

Nur für Forschungszwecke

CD24 Monoklonaler Antikörper

Katalog-Nr.:CL488-65121



Allgemeine Informationen

Katalog-Nr.: CL488-65121	GenBank-Zugangsnummer: BC075622	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 µl/test	GeneID (NCBI): 12484	CloneNo.: M1/69
Wirt: Ratte	Vollständiger Name: CD24a antigen	Anregungs-/Emissionsmaxima-Wellenlängen: 493 nm / 522 nm
Isotyp: IgG2b, kappa		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 mouse splenocytes were surface stained with 5 μ l CoraLite® Plus 647 Anti-Mouse CD3 (17A2) (CL647-65077, Clone: 17A2) and CoraLite® Plus 488-conjugated rat IgG2b isotype control (left) or 5 μ l CoraLite® Plus 488 Anti-Mouse CD24 (CL488-65121, Clone: M1/69) (right). Cells were not fixed.