

Nur für Forschungszwecke

IgG Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL488-65210

Allgemeine Informationen

Katalog-Nr.:
CL488-65210
Größe:
100ug, 0.5 mg/mL
Wirt:
Armenian Hamster
Isotyp:
IgG

GenBank-Zugangsnummer:
GeneID (NCBI):
Vollständiger Name:

Reinigungsmethode:
Affinitätsreinigung
CloneNo.:
PIP
Anregungs-/Emissionsmaxima-Wellenlängen:
493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen

The PIP immunoglobulin is useful as an isotype-matched control. The PIP immunoglobulin is specific for Glutathione S-transferase and can be used as an isotype control for Armenian Hamster IgG antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern. Vor Licht schützen.
Lagerungspuffer:
PBS mit 0,1% Natriumazid und 0,5% BSA, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ BALB/c mouse splenocytes were surface stained with CoraLite® Plus 488-conjugated Armenian Hamster IgG Isotype Control (CL488-65210, Clone: PIP) (purple), or 0.5 ug FITC Anti-Mouse CD3ε (FITC-65060, Clone: 145-2C11) (green), or 0.5 ug FITC-Armenian Hamster IgG Isotype Control (black). Cells were not fixed.