

Nur für Forschungszwecke

# IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL488-65244

## Allgemeine Informationen

Katalog-Nr.:

CL488-65244

Größe:

100ug , 200 µg/ml

Wirt:

Maus

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

MOPC-173

Anregungs-/Emissionsmaxima-  
Wellenlängen:

493 nm / 522 nm

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human PBMCs were surface stained with 0.2 ug Coralite® Plus 488 Anti-Mouse IgG2a Isotype Control (CL488-65244, Clone: MOPC-173) (red), CL488 Anti-Human CD3 (green), or unstained (blue). Cells were not fixed. Lymphocytes were gated.