

Nur für Forschungszwecke

TBR1 Monoklonaler Antikörper

Katalog-Nr.:CL488-66564



Allgemeine Informationen

Katalog-Nr.: CL488-66564	GenBank-Zugangsnummer: BC104844	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 10716	CloneNo.: 5C3A5
Wirt: Maus	Vollständiger Name: T-box, brain, 1	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Isotyp: IgG1	Berechnete Masse: 682 aa, 74 kDa	
Immunogen Katalognummer: AG15262	Beobachtete Masse: 74 kDa	

Anwendungen

Geprüfte Anwendungen:
Getestete Reaktivität:
Hausschwein, Human, Maus, Ratte

Hintergrundinformationen

TBR1, also named as T-box brain protein 1, is a 682 amino acid protein, which contains one T-box DNA-binding domain and localizes in the nucleus. TBR1 is expressed in the brain and as a transcriptional regulator is involved in developmental processes. TBR1 is required for normal brain development. The antibody is conjugated with CL488, Ex/Em 488 nm/515 nm.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

