

Nur für Forschungszwecke

XIAP Monoklonaler Antikörper

Katalog-Nr.:CL488-66800

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:

CL488-66800

Größe:

100ul , Konzentration: 1000 µg/ml von331

Nanodrop;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG18088

GenBank-Zugangsnummer:

BC032729

GeneID (NCBI):

1000

Vollständiger Name:

X-linked inhibitor of apoptosis

Berechnete Masse:

60 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

2A1A12

Empfohlene Verdünnungen:

IF 1:50-1:500

Anregungs-/Emissionsmaxima-

Wellenlängen:

493 nm / 522 nm

Anwendungen

Geprüfte Anwendungen:

IF

Getestete Reaktivität:

Human, Maus, Ratte

Positivkontrollen:

IF : humanes Kolonkarzinomgewebe,

Hintergrundinformationen

XIAP, also named as API3, BIRC4 and IAP3, belongs to the IAP family. It has E3 ubiquitin-protein ligase activity. It mediates the proteasomal degradation of target proteins, such as caspase-3, SMAC or AIFM1. XIAP is an inhibitor of caspase-3, -7 and -9. It mediates activation of MAP3K7/TAK1, leading to the activation of NF-kappa-B. XIAP is an apoptotic suppressor. It is ubiquitinated and degraded by the proteasome in apoptotic cells.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Vor Licht schützen.

Lagerungspuffer:

BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed human colon cancer tissue using Coralite® Plus 488 XIAP antibody (CL488-66800, Clone: 2A1A12) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human colon cancer tissue using Coralite® Plus 488 XIAP antibody (CL488-66800, Clone: 2A1A12) at dilution of 1:200.