

Nur für Forschungszwecke

CD300A Monoklonaler Antikörper

Katalog-Nr.:CL488-67242



Allgemeine Informationen

Katalog-Nr.: CL488-67242	GenBank-Zugangsnummer: BC032352	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von11314 Nanodrop;	GeneID (NCBI): 11314	CloneNo.: 2F9C5
Wirt: Maus	Vollständiger Name: CD300a molecule	Anregungs-/Emissionsmaxima-Wellenlängen: 493 nm / 522 nm
Isotyp: IgG1	Berechnete Masse: 299 aa, 33 kDa	
Immunogen Katalognummer: AG5776	Beobachtete Masse: 50-60 kDa	

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ THP-1 cells were surface stained with 0.4 ug CoraLite® Plus 488 Anti-Human CD300A (CL488-67242, Clone: 2F9C5) (red) or 0.4 ug isotype control. Cells were not fixed.