

Nur für Forschungszwecke

BTN3A3 Monoklonaler Antikörper

Katalog-Nr.:CL488-67560



Allgemeine Informationen

Katalog-Nr.: CL488-67560	GenBank-Zugangsnummer: BC015815	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von10384 Nanodrop;	GeneID (NCBI): von10384	CloneNo.: 1A3B5
Wirt: Maus	Vollständiger Name: butyrophilin, subfamily 3, member A3IF 1:50-1:500	Empfohlene Verdünnungen:
Isotyp: IgG1	Berechnete Masse: 584 aa, 65 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 493 nm / 522 nm
Immunogen Katalognummer: AG12755	Beobachtete Masse: 65 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : Maus-Kolongewebe,
Getestete Reaktivität: Human, Maus	

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glyzerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed mouse colon tissue using CoraLite® Plus 488 BTN3A3 antibody (CL488-67560, Clone: 1A3B5) at dilution of 1:200.