

Nur für Forschungszwecke

METTL3 Monoklonaler Antikörper

Katalog-Nr.:CL488-67733

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.: CL488-67733	GenBank-Zugangsnummer: BC001650	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von56339	GeneID (NCBI): von56339	CloneNo.: 2D8H1
Nanodrop;	Vollständiger Name: methyltransferase like 3	Anregungs-/Emissionsmaxima-Wellenlängen: 493 nm / 522 nm
Wirt: Maus	Berechnete Masse: 64 kDa	
Isotyp: IgG1	Beobachtete Masse: 70 kDa	
Immunogen Katalognummer: AG7110		

Anwendungen

Geprüfte Anwendungen:
Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

METTL3 is a key S-adenosyl-L-methionine-binding subunit, which is component of a complex multicomponent enzyme that catalyzes the methylation of internal adenosine residues in eukaryotic mRNA, forming N6-methyladenosine. It contains 2 putative nuclear localization signals and 2 consensus methylation motifs. The calculated molecular weight of METTL3 is 64 kDa, but modified METTL3 is about 65-70 kDa.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

