

Nur für Forschungszwecke

c-MYC Polyklonaler Antikörper

Katalog-Nr.:CL555-10828

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.: CL555-10828	GenBank-Zugangsnummer: BC000141	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul , Konzentration: 1000 µg/ml von4609	GeneID (NCBI): v-myc myelocytomatosis viral oncogene homolog (avian)	Empfohlene Verdünnungen: IF 1:50-1:500
Nanodrop;	Vollständiger Name: v-myc myelocytomatosis viral oncogene homolog (avian)	Anregungs-/Emissionsmaxima-Wellenlängen: 557 nm / 570nm
Wirt: Kaninchen	Berechnete Masse: 49 kDa	
Isotyp: IgG	Beobachtete Masse: 62-65 kDa, 50 kDa	
Immunogen Katalognummer: AG1263		

Anwendungen

Geprüfte Anwendungen: FC (Intra), IF	Positivkontrollen: IF : HEK-293-Zellen,
Getestete Reaktivität: Human, Maus	

Hintergrundinformationen

MYC contains one basic helix-loop-helix (bHLH) domain. This protein is a multifunctional, nuclear phosphoprotein that plays a role in cell cycle progression, apoptosis and cellular transformation. It functions as a transcription factor that regulates transcription of specific target genes. It seems to activate the transcription of growth-related genes. It binds DNA in a non-specific manner, yet also specifically recognizes the core sequence 5'-CAC[GA]TG-3'. Mutations, overexpression, rearrangement and translocation of this gene have been associated with a variety of hematopoietic tumors, leukemias and lymphomas, including Burkitt lymphoma. This antibody is a rabbit polyclonal antibody raised against recombinant protein of human MYC. The 50kDa band recognized by antibody is the native form of MYC, while the other bands, between 60-70kDa, are the phosphorylated form of MYC [PMID: 12189186]

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using CoraLite®555 c-MYC antibody (CL555-10828) at dilution of 1:200.



1X10⁶ Raji cells were intracellularly stained with 0.2 ug CoraLite®555 Anti-Human c-MYC (CL555-10828) (red), or 0.2 ug Control Antibody. Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).