

Nur für Forschungszwecke

# ITK Polyklonaler Antikörper

Katalog-Nr.:CL555-21525



## Allgemeine Informationen

<b>Katalog-Nr.:</b> CL555-21525	<b>GenBank-Zugangsnummer:</b> BC109077	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 100ul, Konzentration: 1000 µg/ml von3702	<b>GeneID (NCBI):</b> 3702	<b>Empfohlene Verdünnungen:</b> IF 1:50-1:500
<b>Nanodrop;</b>	<b>Vollständiger Name:</b> IL2-inducible T-cell kinase	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 557 nm / 570nm
<b>Wirt:</b> Kaninchen	<b>Berechnete Masse:</b> 620 aa, 72 kDa	
<b>Isotyp:</b> IgG	<b>Beobachtete Masse:</b> 68-72 kDa	
<b>Immunogen Katalognummer:</b> AG16132		

## Anwendungen

<b>Geprüfte Anwendungen:</b> IF	<b>Positivkontrollen:</b> IF : Jurkat-Zellen,
<b>Getestete Reaktivität:</b> Human, Maus, Ratte	

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Vor Licht schützen.  
**Lagerungspuffer:**  
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed Jurkat cells using CoraLite®555 ITK antibody (CL555-21525) at dilution of 1:200.