

Nur für Forschungszwecke

# Cytokeratin 5 Polyklonaler Antikörper



Katalog-Nr.:CL555-28506

## Allgemeine Informationen

<b>Katalog-Nr.:</b> CL555-28506	<b>GenBank-Zugangsnummer:</b> BC024292	<b>Reinigungsmethode:</b> Antigen-Affinitätsreinigung
<b>Größe:</b> 100ul , Konzentration: 1000 µg/ml von3852	<b>GeneID (NCBI):</b> 3852	<b>Empfohlene Verdünnungen:</b> IF 1:200-1:800
<b>Nanodrop;</b>	<b>Vollständiger Name:</b> keratin 5	<b>Anregungs-/Emissionsmaxima-Wellenlängen:</b> 557 nm / 570nm
<b>Wirt:</b> Kaninchen	<b>Berechnete Masse:</b> 590 aa, 62 kDa	
<b>Isotyp:</b> IgG		
<b>Immunogen Katalognummer:</b> AG25436		

## Anwendungen

<b>Geprüfte Anwendungen:</b> IF	<b>Positivkontrollen:</b> IF : A431-Zellen, humanes Tonsillitisgewebe
<b>Getestete Reaktivität:</b> Human, Maus, Ratte	

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Vor Licht schützen.  
**Lagerungspuffer:**  
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Ethanol) fixed A431 cells using Coralite®555 Cytokeratin 5 antibody (CL555-28506) at dilution of 1:400.