

Nur für Forschungszwecke

Beta Catenin Polyklonaler Antikörper



Katalog-Nr.:CL555-51067

Vorgestelltes Produkt

Allgemeine Informationen

Katalog-Nr.: CL555-51067	GenBank-Zugangsnummer: BC058926	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul , Konzentration: 1000 µg/ml von1499	GeneID (NCBI): von1499	Empfohlene Verdünnungen: IF 1:50-1:500
Nanodrop;	Vollständiger Name: catenin (cadherin-associated protein),	Anregungs-/Emissionsmaxima-Wellenlängen: 557 nm / 570nm
Wirt: Kaninchen	beta 1, 88kDa	
Isotyp: IgG	Berechnete Masse: 781 aa, 86 kDa	
	Beobachtete Masse: 92 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : MCF-7-Zellen, humanes Leberkarzinomgewebe
Getestete Reaktivität: Hausschwein, Human, Maus, Ratte	

Hintergrundinformationen

β -Catenin, also known as CTNNB1, is an evolutionarily conserved, multifunctional intracellular protein. β -Catenin is a 92-kDa protein, originally identified in cell adherens junctions (AJs) where it functions to bridge the cytoplasmic domain of cadherins to α -catenin and the actin cytoskeleton. Besides its essential role in the AJs, β -catenin is also a key downstream component of the canonical Wnt pathway that plays diverse and critical roles in embryonic development and adult tissue homeostasis. The Wnt/ β -catenin pathway is also involved in the activation of other intracellular messengers such as calcium fluxes, JNK, and SRC kinases. Deregulation of β -catenin activity is associated with multiple diseases including cancers. (PMID: 22617422; 18334222)

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.

Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Methanol) fixed MCF-7 cells using CoraLite®555 Beta Catenin antibody (CL555-51067) at dilution of 1:100.



Immunofluorescent analysis of (4% PFA) fixed human liver cancer tissue using CoraLite®555 Beta Catenin antibody (CL555-51067) at dilution of 1:100, CoraLite® Plus 647 CD3 antibody (CL647-60181, Clone: 3F3A1, Magenta), CoraLite® Plus 488 CD163 antibody (CL488-16646, green).