

Nur für Forschungszwecke

# IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL555-65208

## Allgemeine Informationen

**Katalog-Nr.:**  
CL555-65208  
**Größe:**  
100ug, 0.5 mg/ml  
**Wirt:**  
Maus  
**Isotyp:**  
IgG2a, kappa

**GenBank-Zugangsnummer:**  
**GeneID (NCBI):**  
**Vollständiger Name:**

**Reinigungsmethode:**  
Affinitätsreinigung  
**CloneNo.:**  
C1.18.4  
**Anregungs-/Emissionsmaxima-Wellenlängen:**  
554 nm / 570 nm

## Anwendungen

**Geprüfte Anwendungen:**  
FC  
**Getestete Reaktivität:**

## Hintergrundinformationen

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

## Lagerung

**Lagerungsbedingungen:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
**Lagerungspuffer:**  
PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



$1 \times 10^6$  mouse splenocytes were surface stained with 0.5  $\mu$ g CoraLite® Plus 555 Mouse IgG2a Isotype Control (CL555-65208, Clone: C1.18.4) (red), or 0.5  $\mu$ g Control Antibody (green), or unstained (blue). Cells were not fixed.