

Nur für Forschungszwecke

# Plvap Monoklonaler Antikörper

Katalog-Nr.: CL555-65214



## Allgemeine Informationen

Katalog-Nr.:  
CL555-65214  
Größe:  
100ug , 500 µg/ml  
Wirt:  
Ratte  
Isotyp:  
IgG2a, kappa

GenBank-Zugangsnummer:  
BC013517  
GeneID (NCBI):  
84094  
Vollständiger Name:  
plasmalemma vesicle associated  
protein

Reinigungsmethode:  
Protein-A-Reinigung  
CloneNo.:  
MECA-32  
Anregungs-/Emissionsmaxima-  
Wellenlängen:  
554 nm / 570 nm

## Anwendungen

Geprüfte Anwendungen:  
FC  
Getestete Reaktivität:  
Maus

## Hintergrundinformationen


### Lagerung

Lagerungsbedingungen:  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.  
Lagerungspuffer:  
PBS with 0.1% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> bEnd.3 cells were surface stained with 1 ug CoraLite® Plus 555 Anti-Mouse Pivap (CL555-65214, Clone:MECA-32) (red), or 1 ug Isotype Control. Cells were not fixed.