

Nur für Forschungszwecke

CD56 Monoklonaler Antikörper

Katalog-Nr.:CL555-65264



Allgemeine Informationen

Katalog-Nr.: CL555-65264	GenBank-Zugangsnummer: BC014205	Reinigungsmethode: Affinitätsreinigung
Größe: 100tests , 5 ul/test	GeneID (NCBI): 4684	CloneNo.: B-A19
Wirt: Maus	Vollständiger Name: neural cell adhesion molecule 1	Anregungs-/Emissionsmaxima-Wellenlängen: 554 nm / 570 nm
Isotyp: IgG1		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:
Human

Hintergrundinformationen

Neural cell adhesion molecule 1 (NCAM1, also known as CD56) is a cell adhesion glycoprotein of the immunoglobulin (Ig) superfamily. It is a multifunction protein involved in synaptic plasticity, neurodevelopment, and neurogenesis. NCAM1 is expressed on human neurons, glial cells, skeletal muscle cells, NK cells and a subset of T cells, and the expression is observed in a wide variety of human tumors, including myeloma, myeloid leukemia, neuroendocrine tumors, Wilms' tumor, neuroblastoma, and NK/T cell lymphomas.

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.
Lagerungspuffer:
PBS with 0.09% sodium azide and 0.5% BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface co-stained with CL750 Anti-Human CD3 and 5 ul CoraLite® Plus 555 Anti-Human CD56 (CL555-65264, Clone: B-A19) or Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.