

Nur für Forschungszwecke

VPS53 Monoklonaler Antikörper

Katalog-Nr.:CL555-67610



Allgemeine Informationen

Katalog-Nr.: CL555-67610	GenBank-Zugangsnummer: BC029560	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von55275 Nanodrop;	GeneID (NCBI): von55275	CloneNo.: 2B12C3
Wirt: Maus	Vollständiger Name: vacuolar protein sorting 53 homolog (S. cerevisiae)	Empfohlene Verdünnungen: IF 1:200-1:800
Isotyp: IgG1	Berechnete Masse: 832 aa, 94 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 557 nm / 570nm
Immunogen Katalognummer: AG29875	Beobachtete Masse: 94 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : HEK-293-Zellen,
Getestete Reaktivität: Human, Maus, Ratte	

Hintergrundinformationen

VPS53 is a component of the Golgi-associated retrograde protein (GARP) complex, also called VFT (VPS fifty-three) complex, composed of VPS51, VPS52, VPS53 and VPS54. GARP complex functions in traffic from endosomes to the trans-Golgi network (PMID: 15878329). GARP proteins interact with RAB proteins and SNARE proteins.

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using CoraLite®555 VPS53 antibody (CL555-67610, Clone: 2B12C3) at dilution of 1:400.