

Nur für Forschungszwecke

ICAM-1 Polyklonaler Antikörper

Katalog-Nr.:CL594-10831

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.: CL594-10831	GenBank-Zugangsnummer: BC015969	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 100ul , Konzentration: 1000 µg/ml von3383	GeneID (NCBI): 3383	Empfohlene Verdünnungen: IF 1:50-1:500
Nanodrop:	Vollständiger Name: intercellular adhesion molecule 1	Anregungs-/Emissionsmaxima-Wellenlängen: 588 nm / 604 nm
Wirt: Kaninchen	Berechnete Masse: 90 kDa	
Isotyp: IgG	Beobachtete Masse: 85-95 kDa	
Immunogen Katalognummer: AG1188		

Anwendungen

Geprüfte Anwendungen: FC (Intra), IF	Positivkontrollen: IF : humanes Lungenkarzinomgewebe,
Getestete Reaktivität: Human	

Hintergrundinformationen

ICAM-1 (CD54) is a 90-kDa transmembrane glycoprotein of the immunoglobulin superfamily and is critical for the firm attachment and transmigration of leukocytes out of blood vessels and into tissues (PMID: 19307690). ICAM-1 is expressed by several cell types, typically on endothelial cells and cells of the immune system, and its expression can be up-regulated by various stimuli, including TNF- α , INF- γ , IL-1 and thrombin (PMID: 3086451; 9694714; 15979056). It is a ligand for LFA-1 and Mac-1, serves as a receptor for rhinovirus, and is one of several receptors used by Plasmodium falciparum (PMID: 2566624; 2538244; 2475784).

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed human lung cancer tissue using CoraLite®594 ICAM-1 antibody (CL594-10831) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human lung cancer tissue using CoraLite®594 ICAM-1 antibody (CL594-10831) at dilution of 1:200.



1X10⁶ Raji cells were intracellularly stained with 0.5 ug CoraLite®594 Anti-Human ICAM-1 (CL594-10831) (red), or 0.5 ug Control Antibody. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).