

Nur für Forschungszwecke

LPAM-1 Monoklonaler Antikörper

Katalog-Nr.:CL594-65101



Allgemeine Informationen

Katalog-Nr.:

CL594-65101

Größe:

100ug, 0.5 mg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC011184

GeneID (NCBI):

16421

Vollständiger Name:

integrin beta 7

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

DATK32

Anregungs-/Emissionsmaxima-
Wellenlängen:

588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse bone marrow cells were surface stained with 0.5 ug CoraLite® 594 Anti-Mouse LPAM-1 (CL594-65101, Clone:DATK32) or 0.5 ug Isotype Control. Cells were not fixed.