

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.:CL594-65208

Allgemeine Informationen

Katalog-Nr.:

CL594-65208

Größe:

100ug , 0.5 mg/ml

Wirt:

Maus

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

GeneID (NCBI):

Vollständiger Name:

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

C1.18.4

Anregungs-/Emissionsmaxima-
Wellenlängen:

588 nm / 604 nm

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Hintergrundinformationen

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 mouse splenocytes were surface stained with 0.5 μ g CoraLite® 594 Mouse IgG2a Isotype Control (CL594-65208, Clone: C1.18.4) (red), or 0.5 μ g Control Antibody (green), or unstained (blue). Cells were not fixed.