

Nur für Forschungszwecke

CD43 Monoklonaler Antikörper

Katalog-Nr.:CL594-66224



Allgemeine Informationen

Katalog-Nr.: CL594-66224	GenBank-Zugangsnummer: BC012350	Reinigungsmethode: Protein-A-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von6693 Nanodrop;	GeneID (NCBI): 2A11D6	CloneNo.: 2A11D6
Wirt: Maus	Vollständiger Name: sialophorin	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG2b	Berechnete Masse: 400 aa, 43 kDa	Anregungs-/Emissionsmaxima-Wellenlängen: 588 nm / 604 nm
Immunogen Katalognummer: AG23595	Beobachtete Masse: 115 kDa	

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : humanes Tonsillitisgewebe,
Getestete Reaktivität: Human	

Hintergrundinformationen

CD43, also named leukosialin, sialophorin, galactoglycoprotein, leukocyte sialoglycoprotein, SPN, and is a mucin-like type I transmembrane protein. In humans, it is expressed in haematopoietic cells, including T lymphocytes, monocytes, granulocytes, natural killer cells, platelets, and haematopoietic stem cells, with the exception of mature erythrocytes and B cell subpopulations. The unglycosylated form of CD43 migrates on SDS-PAGE with a molecular weight of 54 kDa (PMID: 2943740), whereas CD43 from different haematopoietic cell lines displays increased molecular weights since it is glycosylated with O-linked chains that differ in core structure and sialylation.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 CD43 antibody (CL594-66224, Clone: 2A11D6) at dilution of 1:200.