

Nur für Forschungszwecke

HE4 Monoklonaler Antikörper

Katalog-Nr.:CL594-66557



Allgemeine Informationen

Katalog-Nr.: CL594-66557	GenBank-Zugangsnummer: BC046106	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von10406 Nanodrop;	GeneID (NCBI): 10406	CloneNo.: 5B2D9
Wirt: Maus	Vollständiger Name: WAP four-disulfide core domain 2	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG1	Berechnete Masse: 13 kDa	Anregungs-/Emissionsmaxima-Wellenlängen: 588 nm / 604 nm
Immunogen Katalognummer: AG5936		

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : SKOV-3-Zellen,
Getestete Reaktivität: Human	

Hintergrundinformationen

WFDC2, also named as HE4, is a member of the WFDC domain family. The WFDC domain, or WAP Signature motif, contains eight cysteines forming four disulfide bonds at the core of the protein, and functions as a protease inhibitor in many family members. HE4 is expressed in a number of normal tissues, and it is also highly expressed in a number of tumors cells lines, such ovarian, colon, breast, lung and renal cells lines. In some study, HE4 was referred to as a biomarker for ovarian carcinoma.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed SKOV-3 cells using CoraLite® 594 HE4 antibody (CL594-66557, Clone: 5B2D9) at dilution of 1:200.