

Nur für Forschungszwecke

TBR1 Monoklonaler Antikörper

Katalog-Nr.:CL594-66564

1 Publikationen



Allgemeine Informationen

Katalog-Nr.: CL594-66564	GenBank-Zugangsnummer: BC104844	Reinigungsmethode: Protein-G-Reinigung
Größe: 100ul , Konzentration: 1000 µg/ml von10716 Nanodrop;	GeneID (NCBI): 10716	CloneNo.: 5C3A5
Wirt: Maus	Vollständiger Name: T-box, brain, 1	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG1	Berechnete Masse: 682 aa, 74 kDa	Anregungs-/Emissionsmaxima- Wellenlängen: 588 nm / 604 nm
Immunogen Katalognummer: AG15262		

Anwendungen

Geprüfte Anwendungen: IF	Positivkontrollen: IF : Maushirngewebe,
In Publikationen genannte Anwendungen: IF	
Getestete Reaktivität: Hausschwein, Human, Maus, Ratte	
Zitierte Arten: Human	

Hintergrundinformationen

TBR1, also named as T-box brain protein 1, is a 682 amino acid protein, which contains one T-box DNA-binding domain and localizes in the nucleus. TBR1 is expressed in the brain and as a transcriptional regulator is involved in developmental processes. TBR1 is required for normal brain development. The antibody is conjugated with CL594, Ex/Em 593 nm/614 nm.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Aurélie de Thonel	36385105	Nat Commun	IF

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Vor Licht schützen. Nach dem Versand ein Jahr stabil.
Lagerungspuffer:
BS mit 50% Glycerin, 0,05% Proclin300, 0,5% BSA, pH 7,3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CL594-66564 (TBR1 antibody) at dilution of 1:50.



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using CL594-66564 (TBR1 antibody) at dilution of 1:50.